

## Staining retinal cryostat sections

\*PBS instead of Sorenson's can be used

### Acronyms and shorthands:

PFA	Paraformaldehyde
PBS	Phosphate-buffered saline
Sor	Sorenson's buffer
NDS	Normal Donkey Serun
RT	Room temperature
O/N	Over night

### Immunostaining

#### Day1:

1. Thaw slides for **30min** at **RT**
2. Have chamber with wet tissue paper ready
3. Draw a barrier with hydrophobic pen
4. Remove OCT/NEG50 by short incubation (**5min**) in PBS or Sorenson's (use slide mailer if many slides or chamber if only a few slides) at **RT**.
5. Incubate in **5% NDS with 0.5% Triton-X in PBS** (or Sorenson's) for **1h-2h** at **RT** in chamber.
6. Incubate in primary antibody in **5% NDS with 0.5% Triton-X in PBS** (or Sorenson's) **overnight** at **4°C** in chamber.

#### Day2:

7. Wash at least **5x 10min in PBS** (or Sorenson's), shaking in slide mailer
8. Incubate in secondary antibody for **1.5h** at **RT** in **5% NDS with 0.02% Triton-X in PBS** (or Sorenson's) in chamber. Cover with foil to protect from light from this point on.
9. **5x 10' wash in PBS** (or Sorenson's), shaking in slide mailer
10. Mount with DAPI Fluoromount