

## WISH

### Tissue Collection:

1. Dissect embryos in cold RNase free PBS.
2. Fix in 4% PFA in PBS for two hours at 4C.
3. Wash 3x5 in RNase free PBST at RT.
4. Dehydrate for 5min through 25%, 50%, 75% in MeOH/PBST, and x2 in 100% MeOH at RT.
5. Store in 100% MeOH at -20C.

### Day1

1. Re-hydrate through 75%, 50%, and 25% MeOH for 5min each at RT.
2. Wash 3x5 min in PBST at RT.
3. Bleach embryos for 1hr in fresh 6% H<sub>2</sub>O<sub>2</sub> in PBST.
4. Wash 3x5 min in PBST at RT.
5. Incubate in 10µg/ml proteinase K in PBST at RT as follows:

E6.5	4 min
E7.5	5 min
E8.5	6 min
E9.5	7 min
E10.5	8-10 min

6. Wash in 2mg/ml glycine in PBST at RT for 5 min.
7. Wash 3x5 min in PBST at RT.
8. Fix in 4% PFA/0.2% glutaraldehyde in PBST for 20 min.
9. Wash 2x5 min in PBST at RT.
10. Separate embryos in to appropriate wells of a 12 well plate.
11. Wash in 1:1 PBST/hybe solution for 10 min.
12. Replace wash with RT hybe solution, place plate into a 70C oven, allow embryos to warm to 70C (~15 min), and incubate for 2 hrs at 70C.
13. Replace pre-hybe with 3ml hybe/probe at an absolute probe amount of ~75-150ng and incubate at 70C O/N while shaking gently.

### Day2

1. Rinse embryos at 70C in pre-warmed Wash Solution I for 2 min.
2. Wash 2x30 min at 70C in pre-warmed Wash Solution I
3. Wash for 30 min at 65C in pre-warmed Wash Solution I
4. Wash 3x30 min at 65C in pre-warmed Wash Solution II
5. Let embryos cool to RT (about ~20 min)
6. Wash in 3ml of MABT for 3x5 min
7. Block embryos in 10% HISS/10%BBR (Roche 1096176) in MABT at RT for at least 90 min.
8. Pre-absorb antibody:

a) Add 0.5mg of embryo powder to 1ml 2%BBR/MABT and incubate at 70C for 30 min, vortexing briefly every five minutes.

- b) Vortex for 3 min and cool on ice for 2 min.
- c) To the embryo powder/MABT mix add 5 $\mu$ l of HISS and 1 $\mu$ l of anti-DIG AP-conjugate for **each well** of embryos you will have and nutate for at least 1hr at 4C  
(e.g. 5 wells = 25 $\mu$ l HISS and 5 $\mu$ l Ab).
- d) Quick spin to pellet embryo powder and remove the supernatant

#### 9. Prepare antibody mix:

- a) Add 2%BRR/MABT to a 15ml or 50ml tube so you have enough to fill each well to 2.5ml (e.g. 5 wells = 12.5ml).
- b) To this add 20 $\mu$ l of HISS for every well (e.g. 5 wells = 100 $\mu$ l)
- c) Finally, add the full volume of pre-absorbed antibody. This will result in a final [HISS] of 1% and an [AB] of 1:2500.

10. Replace blocking solution with 2.5ml/well of antibody mix and incubate O/N at 4C on a nutator.

### Day3

1. Wash embryos 3x5 with 3ml of MABT at RT
2. Wash all day in MABT at RT changing the MABT every hour or so
3. Wash O/N with MABT at 4C
4. Do another day of washes if necessary

### Day4

1. Wash 3x10 min in AP buffer/1% Levamasole at RT
2. Incubate embryos in staining solution at RT in the dark until desired signal appears
3. Stop reaction with 2mM EDTA for at least 5 min.
4. Fix in 4% PFA/0.1% glutaraldehyde in PBST for 1hr at RT or O/N at 4C
5. Wash 3x5 min in 2mM EDTA/PBST
6. Clear embryos in 80% glycerol/PBST and store at 4C

### Solutions:

#### Stop Buffer

1% SDS  
20mM EDTA  
20mM Tris pH 7.5  
100mM NaCl

#### PFA (80ml):

4% PFA (10ml 32% Electron Microscopy Sciences paraformaldehyde)  
70ml PBS

#### PFA (40ml):

4% PFA(10ml 16% Electron Microscopy Sciences paraformaldehyde)  
30ml PBS

#### 10X PBS (1L):

80g NaCl  
2g KCl

6.1g Na<sub>2</sub>HPO<sub>4</sub> (dibasic,anhydrous)

2g KH<sub>2</sub>PO<sub>4</sub>

800ml dH<sub>2</sub>O

pH 7.4 via NaOH or HCl

DEPC treat and autoclave

**PBST (1L):**

100ml 10X PBS

10ml 10% Triton X-100

bring to volume via DEPC H<sub>2</sub>O

**1M Triethanolamine pH 8.0 (500ml):**

300ml DEPC dH<sub>2</sub>O

66.5ml triethanolamine (sigma T1377 à stock is 7.5M)

pH 8.0 via HCl

final volume 500ml via DEPC treated dH<sub>2</sub>O

**0.1M Triethanolamine pH 8.0 / 0.25 % acetic anhydride(50ml):**

45ml DEPC dH<sub>2</sub>O

5ml 1M triethanolamine pH 8.0

125µl acetic anhydride

**Hybridization Buffer (500ml)**

250ml formamide

0.1% SDS (2.5ml of 20%)

89ml DEPC H<sub>2</sub>O

50µg/ml Yeast tRNA (0.5ml 50mg/ml)

125ml 20X SSC

50µg/ml Heparin (0.5ml 50mg/ml)

60mM Citric Acid (30ml 1M Citric Acid)

0.1% Tween 20 (0.5ml 100%)

**20X SSC (1L):**

3M NaCl (175.3g)

0.3M C<sub>6</sub>H<sub>5</sub>Na<sub>3</sub>O<sub>7</sub>·2H<sub>2</sub>O (dihydrus sodium citrate) (88.23g)

pH 7.0 via HCl

Final volume 1L autoclaved dH<sub>2</sub>O

DO NOT AUTOCLAVE

**5X MABT (1L):**

58g Maleic Acid

43.5g NaCl

pH 7.5 with Tris Base (need ~100g)

final volume 1L via autoclaved dH<sub>2</sub>O

**Alkaline Phosphatase Buffer (100ml):**

100mM Tris pH 9.5 (10ml of 1M)

50mM MgCl<sub>2</sub> (5ml of 1M)

100mM NaCl (2ml of 5M)

0.1% Tween 20 (1ml of 10%)

**Wash Solution I (500ml):**

250ml formamide  
25ml 20% SDS  
70ml DEPC H<sub>2</sub>O  
125ml 20X SSC  
30ml 1M Citric Acid

**Wash Solution II (500ml):**

250ml formamide  
5ml 20% SDS  
182.5ml DEPC H<sub>2</sub>O  
50ml 20X SSC  
12ml 1M Citric Acid  
0.5ml Tween 20

12 well plates are from EMS (cat# 71549-08)